

micra

Digital Diagnostic Audiometer (with USB)



Features :-

Tests	Air Conductions, Bone Conductions and Speech Audiometry
Operation	Soft Feather Touch Control Keys
Interrupter	Normal and Pulse Tone
Free field Test	Built in FF Amplifier
PC connection	Thru USB Port
Frequency Range	125 Hz To 12000 Hz
Intensity Range	AC - (-) 10dB To 120 dB
	BC - (-) 10dB To 80 dB
	Speech - (-) 10dB To 100 dB
	FF - (-) 10dB To 120 dB
Distortion	Less than 0.1% in Tone
Special Feature	Tone decay, Talk over and Talk back
Masking	Wide Band
Display	16 Characters 2 Row LCD Display
Weight	1500 gms
Calibration	Software Calibration
Power Supply	AC/DC Adaptor 12 V -2Amp DC Output
Standard Accessories	DR 59 Head Phones, B-71 Radio Ear Bone Vibrator with Headband, Talk Back / Talk Forward MIC, Monitor Phone, Patient Response Button, 12 V DC Power Supply, Software CD, Carrying Case, USB Cable and Audiogram Pad

AUDITIVO HEARING SERVICES PVT. LTD.

Flat No. 309, 19, Oriental Apartment, Commercial Centre, Yusuf Sarai, New Delhi- 49
info@auditivohearing.in, sales@auditivohearing.in | +91 9810112723, 9810112042